Vend Test Database

	11/01/200 500ml	1	*Model: *Type:	501E SERIES 90 PET		Test #: 20 Country: *Description	CHILE		TURAL
Pkg. Dia. Pkg. Dia.	. (inches): . (mm):	2.39 60.71		Pkg Ht. (inches): Pkg. Ht. (mm):	9.34 237.24	● New ○ Revi		Instructions	5:

ECN #: 2001-183-001

Wide Co	Wide Column Setup Deep: 2							
Qty	Part Number	Description						
4	801,812,04x.x1	3/32 flat shim (2 each side)						
2	610,070,06x.x3	shim angle (1 each side on top of flat shim next to package)						
Notes:								

Cam set in hole #: 3 Rear Spacer set in hole #: 24 Load Package: caps to the back

Narrow Column Setup Deep: 2

Qty	Part Number	Description
1	805,701,14x.x1	can stop shim - right
1	609,072,80x.x3	bottle neck guide assy.
1	432,070,27x.x3	large product sold out paddle
1	801,903,24x.x1	white tube
1	801,903,23x.x1	red tube
2	800,503,83x.x1	Assy, R, Rod bottle E-Model

Notes: Attach shim angle to top of can stop shim. Place one rod with red tube in hole "A" of rotor. Place one rod with white tube in hole "C" of rotor.

Cam set in hole #: 5Rear Spacer set in hole #: 24Load Package: caps to the rear

Discharge Member:large portWide ColumnClosure Strip:large portNarrow ColumCan Chute Accessories:large portTotal Capacity

Wide Column Capacity:62Narrow Column Capacity:30Total Capacity:302