Vend Test Database

ECN #: 002042

Date: () *Size: { Release Level:)7/02/2004 500 ML	4	*Model: *Type:	501E SERIES 90 PET	1	Coun	·	SINGAP	
Pkg. Dia. Pkg. Dia.	(inches): (mm):	2.6 66.04		Pkg Ht. (inches): Pkg. Ht. (mm):	9.34 237.24		● New ○ Revise	d	Instructions:

Qty	olumn Setup Deep: 2 Part Number	Description		
-	801,812,04x.x1	3/32 flat shim (2 each side next to column wall)		
	801,809,81x.x1	3/32 short flat shim (2 each side next to bottles)		
	801,809,80x.x1	vending cam, grey		
1	801,809,79x.x1	adjustable cam, red		
Notes:				

Cam set in hole #: 3 Rear Spacer set in hole #: 24 Load Package: caps to the back

Narrow Column Setup Deep: 2

Qty	Part Number	Description
1	805,701,15x.x1	bottle shim, right
1	801,903,72x.x1	assembly, bottle divider
1	800,503,83x.x1	rod & spring assy

Notes: Place rod in hole "B" of rotor.

Cam set in hole #: 2 Rear Spacer set in hole #: 24 Load Package: caps to the back

Discharge Member:large portClosure Strip:large portCan Chute Accessories:large port

Wide Column Capacity:56Narrow Column Capacity:30Total Capacity:296