## **Vend Test Database**

ECN#: 001924 Test #: 97-531-08

Date: 05/26/2004 \*Model: 501E SERIES 90 Country: **CANADA** \*Size: 600 ML \*Type: PET \*Description: BARQ'S

Release Level:

Pkg. Dia. (inches): 2.89 Pkg Ht. (inches): 9.05

Pkg. Dia. (mm): Pkg. Ht. (mm): 73.4 229.9 O New Revised Instructions:

Wide Column Setup Deep: 2			
Qty	Part Number	Description	
2	801,809,81x.x1	3/32 short flat shim (1 each side)	

Notes:

Cam set in hole #: 3 Rear Spacer set in hole #: 24

Load Package: caps to the back

Narrow Column Setup Deep: 2

Qty	Part Number	Description
1	805,701,15x.x1	bottle shim, right
1	609,072,80x.x3	bottle neck guide assy.
2	800,902,40x.x1	vertical metal guide clip
1	801,810,57x.x1	front guide strip

Notes: Assemble bottle neck guide assy, to rear spacer. Install vertical guide clips in left column wall. Hang front guide strip on clips.

Cam set in hole #: 5 Rear Spacer set in hole #: 24

Load Package: Back bottles first caps to the back, then front bottles caps to the back.

Discharge Member: Wide Column Capacity: large port Closure Strip: Narrow Column Capacity: 26 large port Can Chute Accessories: large port **Total Capacity:** 258