Vend Test Database

ECN #: 001933 Test #: 97-545-10

 Date:
 05/27/2004
 *Model:
 276E SERIES 90
 Country:
 ENGLAND

 *Size:
 500 ML
 *Type:
 PET
 *Description:
 FANTA

Release Level:

Pkg. Dia. (inches): 2.68 Pkg Ht. (inches): 9.18

Pkg. Dia. (mm): 68 Pkg. Ht. (mm): 233.2

\bigcirc	New
	Revised

Instructions:

Wide Column Setup Deep: n/a		
Qty	Part Number	Description
		·

Notes:

Cam set in hole #: Rear Spacer set in hole #:

Load Package:

Narrow Column Setup Deep: 2

Qty	Part Number	Description
1	805,701,15x.x1	bottle shim, right
1	609,072,80x.x3	bottle neck guide assy.
1	609,010,50x.x3	assy. column wall, right
1	800,503,83x.x1	Assy, R, Rod Bottle E-Model
1	801,903,25x.x1	yellow tubing

Notes: Place rod with yellow tube with front in hole "B" and rear in hole "A" of rotor. Install column wall assy. Assemble bottle neck guide assy. to rear spacer.

Cam set in hole #: 5 Rear Spacer set in hole #: 24

Load Package: caps to the back

Discharge Member:large portWide Column Capacity:Closure Strip:large portNarrow Column Capacity: 32Can Chute Accessories:large portTotal Capacity: 224